


# Extracting DNA from a banana


Authors: Victor Sojo, Aishwaria Lohi | References: [youtu.be/ew9-YGrpWo](https://youtu.be/ew9-YGrpWo) | [bit.ly/2AJmwqI](https://bit.ly/2AJmwqI)

v02



## 1. Chop the banana

Cut half of a banana (ripe is better), and peel it. Then chop it finely. Transfer it to cup #1 where you will mash it. You may eat the other half ☺




## 2. Mash it very well

Use a pestle if you have one, otherwise a spoon will do. Spend a few minutes doing this. We need to make sure the tissues of the banana are broken down well.



## 4. Add the lysis solution and mix very well

Pour the lysis solution onto the banana paste, and mix very well, for at least 5 minutes. You can use the pestle again if you have it, otherwise the spoon again




## 3. Make "lysis" solution

DNA is trapped deep inside the cells of the banana. We need to break it free (this is called "cell lysis", breaking the cell). In cup #2, mix about 1 teaspoon of washing up liquid (or liquid soap) and some 10 tablespoons of water (no need to be too exact here).



## 5. Filter the solution

Using a sieve (sifter/strainer) with a fine mesh, extract the liquid into cup #3. You can also use a coffee filter, or a funnel with a layer or two of toilet or kitchen tissue.




## 6. Transfer to a tall transparent container

The ideal container would be a big test tube, but if you don't have one, try to get the thinnest glass jar you can find.

## 7. Extract the DNA


Lean the container to around 45° and add alcohol at the top, letting it flow down the walls slowly. About a finger in height is enough. Almost immediately, you will see a white slime at the interface. That's DNA! It "precipitates" because the water in which it was dissolved goes with the alcohol instead.



## 8. Fish the DNA out

Using a stick appropriate to the shape of your container, try to fish out the DNA. Note that this is not pure DNA, it has RNA, proteins, and a lot of other things. We would need to purify it with professional reagents if we wanted to analyse it. You can try a similar method with tomatoes, strawberries, and even saliva! (starting at step 3).

We recommend that you now do the *DNA & Genetics* activity.



## Materials, samples and reagents

- Banana (or tomato, strawberries)
- Alcohol (as pure as you can find)
- Washing-up liquid (soap)
- Water
- 3 cups or small glasses
- Knife (doesn't have to be sharp)
- Teaspoon & Tablespoon
- Pestle (optional)
- Sieve/sifter (or coffee filter)
- Test tube (or tall thin glass/jar)
- Saliva (optional)



[science-corps.org](https://science-corps.org)



[aavishkaar-palampur.org](https://aavishkaar-palampur.org)



[eduversum.org](https://eduversum.org)



This content is released under a CC-BY-NC-SA license. In short, this means you can use and modify it freely as long as you: do NOT sell what you create, you share it under the same conditions, and you credit us by linking to our website(s), where you'll find additional information.